



Coag2D

**THE MOST VERSATILE 2 CHANNEL
OPTICAL COAGULOMETER
AT YOUR HAND**

**small compact design
supports D-dimer and AT-III test
ready for coagulation, chromogenic and turbidimetric tests**



Coag2D

COAG 2D COAGULOMETER SPECIFICATIONS

Instrument features:

- 2 independent measuring channels
- 2 reagent positions (2 mL, 5mL)
- 4 incubation positions
- Coagulation screening test
- Quantitative D-dimer measurement (Coag 2D PLUS only)
- AT-III activity measurement (Coag 2D PLUS only)
- QC feature
- Supports application of Master Calibration Curve
- Maintenance free
- Up to 1000 results storage
- Extendable storage capacity
- Connection to Coag Data data management software

System elements:

- 2 thermostated reagent positions with magnetic stirring
- Colour LCD touchscreen
- Innovative user interface
- LIS host interface connection
- Small, compact footprint, easy too operate
- 640 nm wavelength detection
- 405 nm and 570 nm wavelength detection (Coag 2D PLUS only)



EUROPE

Diagon Ltd. - Hungary
1047 Budapest, Baross u. 48-52.
tel.: (+36 1) 369 6500
fax: (+36 1) 369 6301
e-mail: diagon@diagon.com

www.diagon.com

SOUTH AMERICA

DIAGON Ltda. - Brasil
Av. Joaquim José Diniz, 787
CEP 31910-520
Bairro Fernão Dias
Belo Horizonte/MG
Brasil
Tel: (+55 31) 3488 6777
Fax: (+55 31) 3488 6778
e-mail: diagon@diagon.com.br
www.diagon.com.br



ASIA

PT. Inti Diagonama Selaras - Indonesia
Perkantoran Plaza Pasifik A3/59 -
Jl. Raya Boulevard Barat
Kelapa Gading Permal Jakarta 14250
Indonesia
tel.: (+62 21) 458 40649, 458 40650
fax: (+62 21) 458 45024
e-mail: dsu@centrin.net.id